

# DENTAPEN™

## Dental anesthesia: Painless Cordless Effortless



Compact handheld device designed for dentists to provide painless anesthesia

# DENTAPEN™ cutting-edge C-CLAD\* technology for effortless dental anesthesia

Dentapen offers patients better comfort by significantly reducing pain and eliminating fear of injections.

Dentapen helps dentists avoid muscle strain caused by repeated manual injections.

Suitable for all anesthesia techniques: palatal, infiltration, block and intraligamental.

## Two-way activation

Cordless, Battery powered

Compact and light (40g)

Compatible with standard 1.7mL cartridges

Rotating barrel for easy bevel placement

Self-aspiration

Autoclavable cartridge holder

3 speeds  
2 modes

8cm



## Two-way activation

Index or thumb activation for better comfort



☀ Slow 90 s/mL

☀ Medium 60 s/mL

☀ Fast 30 s/mL



☀ Intraligamentary mode


☀ Ramp-up mode = anesthesia flow increases gradually for totally painless procedure

\* computer-controlled local anesthesia device

# Patients feel less pain with

## DENTAPEN™

Controlled initial force and constant flow for smoother injections



LESS

IS

MORE



PAIN



FEAR



CHILDREN  
FRIENDLY



SATISFIED  
PATIENTS



LOYALTY



PRESTIGE

Dentapen Kit ..... Ref. DK01

- Dentapen Device x1
- Finger Grip x5
- Cartridge Holder x5
- Battery set



Protective sleeves available separately

J U V A +  
I N J E C T I O N   S Y S T E M S

Manufacture: Maladière 71C, 2002 Neuchâtel, Switzerland  
Sales: Av. de Miremont 12, 1206 Geneva, Switzerland

info@dentapen.ch +41 22 588 66 89

www.dentapen.ch

